



ANIMAL CARE

PRODUCT SPECIFICATION

REF 510002410 SNÖGG ANIMAL FOAM BAND ADHE 9CMX2M 20X1

1 GENERAL PRODUCT DESCRIPTION

Description

Adhesive Foam Bandage is a latex free and elastic foam bandage with adhesive on one side that sticks to fur. Ideal for padding and fixation and must be used in combination with a compress or self-adherent foam bandages. Use Snøgg Removal Wipes or skin friendly oil to ease adhesive removal.

Storage information

Dry storage in room temperature. For maximum protection, keep in original package until use.



2 DETAILED PRODUCT INFORMATION

Content per consumer unit

Product/Material

Foam Bandage

Quantity

9cmx2m

Dimensions (w x l mm)

Bandage: 9cmx2m

Thickness: 5 mm

Material

Bandage: White polyurethane foam without latex.

Box: Printed cardboard

Shelf life: 5 years from production date

3 LOGISTIC INFORMATION

	CU	SKU	Pallet
GTIN:	7310617310984	7310800310364	7312020317369
Content:	1 Roll	20 CU	6 SKU/layer / 3 layers
Size (w x l x h mm):	125 x 92 x 130	493 x 273 x 275	1200 x 800 x 964
Gross weight:	85 g	2000 g	61 kg
Packing material:	Carton 31 g (0 % recycled fibers)	Corrugated Board 300 g (78 % recycled fibers)	-

Version: 09

Issued: 2023-04-06