



ANIMAL CARE

PRODUCT SPECIFICATION

REF 510002379 SNÖGG ANIMAL FOAM BAND THIN 6CMX4,5M 10X1

1 GENERAL PRODUCT DESCRIPTION

Description

The product is a cohesive polyurethane foam bandage with a polymer emulsion adhesive. The foam is rolled with a protective paper on each side. The product is intended to be used as a mechanical barrier and for absorption of exudates from minor wounds.

Storage information

Dry storage in room temperature. For maximum protection, keep in original package until use.



2 DETAILED PRODUCT INFORMATION

Content per consumer unit

Product/Material

Foam Bandage

Quantity

4,5m

Dimensions (w x l mm)

Bandage: 6cmx4,5m

Thickness: 1,9 mm

Material

Bandage: Beige polyurethane foam without latex.

Box: Printed cardboard

Shelf life: 5 years from production date

3 LOGISTIC INFORMATION

	CU	SKU	Pallet
GTIN:	7310617311028	7310800310401	7312020317406
Content:	1 Roll	10 CU	10 SKU/layer / 8 layers
Size (w x l x h mm):	117 x 34 x 115	334 x 245 x 130	1200 x 800 x 1160
Gross weight:	69 g	840 g	92 kg
Packing material:	Carton 26 g (0 % recycled fibers)	Corrugated Board 150 g (78 % recycled fibers)	-

