



ANIMAL CARE

PRODUCT SPECIFICATION

REF 510002411 SNÖGG ANIMAL FOAM BAND THICK 9CMX2M 20X1

1 GENERAL PRODUCT DESCRIPTION

Description

A cohesive polyurethane foam bandage with a polymer emulsion adhesive. The foam is rolled on a paper core with a protective paper on one side. The paper end of the roll is glued for fixation. The product is intended to be used as a mechanical barrier and for absorption of exudates from minor wounds.

Storage information

Dry storage in room temperature. For maximum protection, keep in original package until use.



2 DETAILED PRODUCT INFORMATION

Content per consumer unit

Product/Material

Foam Bandage

Quantity

9cmx2m

Dimensions (w x l mm)

Bandage: 9cmx2m

Thickness: 5 mm

Material

Bandage: Beige polyurethane foam without latex.

Box: Printed cardboard

Shelf life: 5 years from production date

3 LOGISTIC INFORMATION

	CU	SKU	Pallet
GTIN:	7310617311011	7310800310395	7312020317390
Content:	1 Roll	20 CU	6 SKU/layer / 3 layers
Size (w x l x h mm):	125 x 92 x 130	493 x 273 x 275	1200 x 800 x 964
Gross weight:	72 g	1740 g	56 kg
Packing material:	Carton 30 g (0 % recycled fibers)	Corrugated Board 150 g (78 % recycled fibers)	-

Version: 09

Issued: 2023-04-06