



ANIMAL CARE

PRODUCT SPECIFICATION

REF 510002378 SNÖGG ANIMAL CARE REMOVAL WIPE 24X30 PCS

1 GENERAL PRODUCT DESCRIPTION

Description

Folded paper tissues, which are saturated with an adhesive remover substance (100% disiloxane). In order to prevent leakage, the tissues are packed in a well-sealed foil. The sachets are packed in CU boxes, that are packed in SKU boxes. Intended to be used as an adhesive remover.

Storage information

Dry storage in room temperature. For maximum protection, keep in original package until use.

Highly flammable liquid and vapor. Very toxic to aquatic life with long lasting effects. Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Store in a well-ventilated place. Keep coop. Avoid release to the environments. Collect spillage. Keep out of reach for children.



2 DETAILED PRODUCT INFORMATION

Content per consumer unit

<u>Product/Material</u>	<u>Quantity</u>	<u>Dimensions (w x l mm)</u>
Removal Wipe	30 pcs	60 x 80 (Sachet) 120 x 180 (Tissue unfolded)

Material

Wipe:	Non-woven material
Liquid:	Solution with 100 % Hexamethyldisiloxane
Sachet:	Heat sealed laminate paper
Box:	Printed cardboard

Shelf life: 4 years from production date

3 LOGISTIC INFORMATION

	CU	SKU	Pallet
GTIN:	6417495400743	6417495400750	6417495400767
Content:	30 Sachet	24 CU	12 SKU/layer / 3 layers
Size (w x l x h mm):	160 x 62 x 82	384 x 190 x 336	1200 x 800 x 1152
Gross weight:	183 g	4 851 g	97 kg
Packing material:	Cardboard box: 20 g (0 % recycled fibers)	Corrugated Board 467 g (75 % recycled fibers)	-

Version: 09
Issued: 2023-04-06